

# Eye irritation protocol

## Time-To-Toxicity test method

### Liquid raw materials

#### DESCRIPTION

Model: SkinEthic™ HCE

Receipt: Transfer epithelium from agarose to maintenance medium in 6-well plate and incubate at least overnight (37°C)



Transfer tissues to fresh maintenance medium in the same 6-well plates



Preparation of the 20% dilution (weight/volume) in distilled water of methyl acetate (Positive Control) and test chemicals for time treatment 16 and 120min



Treatment: 2 tissues each with 80 µL PBS without Ca<sup>2+</sup> & Mg<sup>2+</sup> (Negative Control) or 10µL PBS- + 80 µL methyl acetate (Positive Control) or 10µL PBS- + 80 µL test chemical (test treatment)

Neat for 5min time treatment or 20% w/v for 16 and 120min time treatment



Treatment Period: Incubate for 5min ± 15s (RT) or for 16min ± 1min (RT) or for 120min ± 2min (37°C)



Rinse with PBS- (25 mL: 2 mL/jet)



Post-Soak Immersion: Immerse tissues in 4 mL fresh maintenance medium in 12-well plates  
Post-Soak Period: Incubate for 10min (RT)



Viability: Transfer tissues into MTT solution in 24-well plates



Incubate tissues for 3h ± 15min (37°C)



Rinse the MTT solution with PBS-



Extraction: Immerse the inserts in 1.5mL isopropanol (formazan extraction) in 24-well plates




Extract formazan (minimum 2h, RT or overnight, 4°C)



Perforate the insert and homogenize formazan extract



Read OD with microplate spectrophotometer at 570 nm



#### PREDICTION MODEL

<i>In vitro</i> Result			Classification ( <i>In vivo</i> Prediction)
<i>Test chemical neat</i>	<i>Test chemical diluted at 20%</i>		
<i>Time treatment: 5 min</i>	<i>Time treatment: 16 min</i>	<i>Time treatment: 120 min</i>	
Mean tissue viability ≤ 50 %	≤ 50 %	≤ 50 %	Category 1
Any other combination of values			Category 2
> 50%	> 50%	> 50%	No Category